



# Sanjeevani

*Sanskrit: 'that which gives life'*

Neuro-Symbolic AI for Automated  
CKD Triage & EMR Drafting via 1D NMR

Abhiram Radha Krishna · Kaggle MedGemma Impact Challenge · 2025

## The Problem in Numbers

850M

people with CKD

90%

undiagnosed early

3hrs

daily EMR burden

0

LLMs can read NMR

*Built entirely from scratch  
for this challenge.*

Google's HAI-DEF models deployed  
strictly within their competence domain.

# A Silent Epidemic & An Interpretability Bottleneck



## The Spectroscopic Bottleneck

<sup>1</sup>H NMR **spectroscopy** of urine is the gold standard for early CKD metabolic profiling. It detects:

**Citrate** depleted in renal tubular dysfunction

**TMAO** elevated uremic toxin

**Lactate** elevated in renal hypoxia

**Creatinine** normalisation anchor metabolite

*BUT: requires highly trained spectroscopists — a critical bottleneck for under-resourced clinics.*



## Why Pure LLMs Fail Here

LLMs are text-prediction engines. They cannot perform spectral deconvolution on raw float arrays.

✗ Hallucinate diagnoses from random noise

✗ Cannot extract metabolite ratios reliably

✗ No deterministic mathematical grounding

✗ Dangerous false confidence on bad data

**Solution: Use LLMs only for what they're good at. Build a deterministic physics layer first.**

# Three AI Agents — Each Doing Only What It Can Do Well



STAGE 1

**PaliGemma***paligemma-3b-mix-224*

## Visual Quality Control

Renders 2.5–4.5 ppm anchor region as image. Classifies 'sharp peaks' vs 'flat noise'. Halts pipeline on degraded samples.



STAGE 2

**SanjeevaniEngine***Custom PyTorch CNN-Transformer*

## Biomarker Deconvolution

5-channel independent decoders extract creatinine, citrate, lactate, TMAO, taurine. Outputs deterministic risk score.



STAGE 3

**MedGemma***medgemma-1.5-4b-it*

## Clinical Communication

Receives verified ratios only — never raw data. Dual persona: Physician Assessment & Plan + patient-friendly EMR draft.



**Safety Gate:** MedGemma never sees raw spectral data — only three verified decimal ratios from the physics engine. This prevents hallucinated clinical reasoning.

# Patient C: Degraded Sample → Pipeline Halts Safely

## Without Sanjeevani

1. Raw noise fed directly to MedGemma
2. LLM pattern-matches text from noise
3. Generates a plausible-sounding diagnosis
4. Clinician acts on a hallucinated result
5. Patient harmed — false negative discharge

## With Sanjeevani

### ① PaliGemma Vision Gate

Renders anchor region (2.5–4.5 ppm). Classifies image as 'flat noise' — not 'sharp peaks'.

### ② Mathematical Creatinine Gate

Checks creatinine anchor amplitude  $< 0.01$  threshold. Secondary biochemical validation — dual-layer QC.

### ③ Pipeline Halts

MedGemma is never invoked. Dashboard shows: 'INCONCLUSIVE (QC FAILURE) — Biomarker Math Halted: Data integrity compromised.'

[ See dashboard screenshot: [inconclusive\\_fig\\_1.png](#) + [inconclusive\\_fig\\_2.png](#) ]

# Patient B: Type 2 Diabetic → HIGH RISK DETECTED

## Biomarker Panel

Creatinine **0.90** (anchor)

*reference peak — normal*

Citrate Ratio **0.06** / Cr

*⚠ severely depleted — renal tubular dysfunction*

TMAO Ratio **1.14** / Cr

*⚠ elevated — uremic toxin accumulation*

Lactate Ratio **1.36** / Cr

*⚠ elevated — anaerobic metabolism*

**RISK SCORE:**

**3.22** Threshold: 1.2 →  
HIGH RISK



## MedGemma — Physician View

*"The depleted citrate ratio, elevated TMAO ratio, and elevated lactate ratio are indicative of metabolic acidosis and potential kidney dysfunction. These findings suggest impaired renal bicarbonate excretion... A renal ultrasound is recommended."*

## MedGemma — Patient EMR Draft

*"Hi [Name], I'm following up on your recent check-up. Your fatigue, along with recent changes in your kidney function, is something we need to look into further... We're recommending a renal ultrasound."*

# SanjeevaniEngine: Architecture & Performance Evolution



## Architecture

Encoder	2× InceptionBlock (kernels 1,3,5,7) + MaxPool1D → 64-channel latent
Sequence	2-layer TransformerEncoder (d=64, heads=4) for long-range dependencies
Decoder	5 INDEPENDENT branches — one per biomarker, with BatchNorm + weighted loss
Loss weights	TMAO 5×, Lactate 3×, Citrate 2×, Taurine 2×, Creatinine 1×
Training	10,000 synthetic pseudo-Voigt spectra · 25 epochs · Adam optimizer



## Performance Evolution


Version	Acc	Sens	Spec
V1 Baseline (shared decoder)	93%	96.6%	87.8%
V2 Indep. decoders + decoys	91%	96.1%	85.7%
Final: weighted loss (all 5 ch. active)	89.5%	100%	76.9%

↓ Accuracy intentional: reflects removing inflated closed-world metrics & forcing all 5 channels active

[ See: [confusion\\_matrix.png](#) — 100% sensitivity, 14 QC flags / 100 patients ]

# MTBLS1 Zero-Shot Transfer: Honest Domain Gap Analysis

## Domain Gap Analysis Figure



## Two Domain Gaps Identified

### Gap 1: Amplitude Mismatch

Synthetic: creatinine always dominant (0.8–1.2).  
Real urine: creatinine can be subdominant. Model outputs zero for creatinine channel.

### Gap 2: Elevated Baseline Floor

MTBLS1 spectra show macromolecular background at 0.25–0.35 normalised intensity. Synthetic data has near-zero baseline. All real peaks compressed.

### V2 Fixes:

(1) Add broad Lorentzian baseline augmentation to generator · (2) Calibrate creatinine amplitude distribution from real MTBLS1 · (3) Spinach quantum spin Hamiltonian simulator for J-coupled multiplets

# Sanjeevani

Neuro-Symbolic AI for Automated CKD Triage & EMR Drafting via 1D NMR



## Dual-Layer QC

PaliGemma vision + creatinine biochemical gate prevent bad data reaching clinical reasoning



## Channel Collapse Fix

Independent decoder branches + weighted loss forces all 5 biomarker channels to function



## Hallucination Prevention

MedGemma receives only verified decimal ratios — never raw spectral data



## Honest Transfer Eval

Zero-shot MTBLS1 test with quantified domain gap analysis and specific V2 roadmap

*"A model that fails honestly and explicably is more valuable for clinical AI than one that succeeds silently on its own test set."*